Convenience Pack - Terminally sterilized
DFU/ Bruksanvisning Healon GV™ PRO, EN/ES/PT/EL
(English/Spanish/Portuguese/Greek)
(engelska/spanska/portugisiska/grekiska)

Format: 420 x 297 mm (A3-format)
Folding/Falsad: 210 x 49,5 mm

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0344

© Johnson & Johnson Surgical Vision, Inc. 2021

Healon GV™ PRO

Ophthalmic Viscosurgical Device

Product information

Healon GV™ PRO viscoelastic solution is for intraocular use. It is terminally sterilized by steam sterilization in disposable syringes delivering 0,85 ml. The syringe is packaged together with a single use 27G canula in a blister which is terminally sterilized by ethylene oxide.

Composition per ml
Sodium hyaluronate SE 18 mg
Phosphate buffer ad 1 ml

Description

The Healon GV PRO Ophthalmic Viscosurgical Device (OVD) is a steam-sterilized, non-pyrogenic solution of a highly purified, non-inflammatory, high molecular mass fraction of sodium hyaluronate (HA) derived by microbial fermentation. The HA is provided at a concentration of 18 mg/ml (1,8%) dissolved in a physiological buffer. The Healon GV PRO OVD has an osmolality and a pH similar to aqueous humour. HA in the Healon GV PRO OVD has an average molecular mass of 3.2 million. The Healon GV PRO OVD has a viscosity at rest of about 2000 Pas, 2 million times higher than that of aqueous humour. At high shear rates, such as during injection, the viscosity of the Healon GV PRO OVD decreases dramatically due to shear thinning, facilitating injection through a 27G canula. The graph represents the flow curve (shear viscosity versus shear rate).

HA is a polysaccharide made up of repeating disaccharide units linked by glycosidic bonds. HA occurs naturally on the corneal endothelium, bound to specific receptors. The natural HA layer on the corneal endothelium is reduced during irrigation, but can be restored by exogenous HA. The Healon GV PRO OVD acts as a scavenger by neutralizing the free radicals formed during cataract surgery using ultrasound.

The Healon GV PRO OVD is easy to inject as the viscosity decreases considerably when expelled through a canula. The Healon GV PRO OVD regains its high viscosity immediately after injection. The Healon GV PRO OVD creates and maintains a deeper anterior chamber than the Healon™ PRO OVD, allowing safe and controlled manipulation, reducing the risk of trauma to the corneal endothelium and other tissues. The Healon GV PRO OVD is used to atraumatically manoeuvre, control and separate tissues.

The Healon GV PRO OVD protects against mechanical impacts. The Healon GV PRO OVD provides counter pressure in cases with high pressure on the globe and can be used to push back an iris prolapse or a bulging vitreous. The Healon GV PRO OVD cretes and maintains a clear field of vision during surgery.

The Healon GV PRO OVD is eliminated in the same way as endogenous HA in the anterior segment, i.e. via the trabecular meshwork. Upon entering the general circulation HA is metabolized in the liver.

Indication

Eye surgery.

Contraindications

None known at present.

Dosage and methods of administration

The Healon GV PRO OVD is injected through a fine canula into the eye. The amount used varies with the type of surgery, repeated injections may be required. The total amount introduced is usually less than 1 ml in anterior segment surgery.

Undesirable effects

A transient rise in intraocular pressure may be seen postoperatively due to a mechanical obstruction of the trabecular meshwork. The peak of this pressure rise may be expected 4-8 hours postoperatively. To minimize this risk it is recommended that the Healon GV PRO OVD is irrigated/aspirated out of the eye at the close of surgery. If complete removal is not achieved the risk of increased intraocular pressure may be higher with the Healon GV PRO OVD than with Healon PRO.

Incompatibility

Sodium hyaluronate solution, like Healon GV PRO OVD, may appear cloudy or form precipitates when it is injected. The physician should be aware of this phenomenon and, should it be observed, remove the cloudy or precipitated material by irrigation and/or aspiration. Based on *in-vitro* laboratory studies, this phenomenon may be related to interactions with concomitantly used ophthalmic medications or detergents which remain in reused canulas. *In-vitro* studies have also shown incompatibility, resulting in opalescence, between sodium hyaluronate and solutions containing cationic components, e.g., detergents, quaternary ammonium compounds and benzalkonium chloride.

Precautions

For single use only. Do not resterilize syringe or canula. Reuse of product or canula may cause inflammatory reaction. Do not use if packaging has been opened or damaged. Use only the canula provided in the package. Use only if solution is clear.

Do not use after date shown on product. Store cold at 2 to 8°C. Protect from light.

Precautions normally considered during ophthalmic surgical procedures should be taken.

Before initiating phacoemulsification, use irrigation/ aspiration to create a fluid filled space above the lens. This reduces the risk of initial visco-occlusion of the phaco tip or the irrigation line which could cause phaco tip heating.

Postoperative intraocular pressure (IOP) may be increased if the Healon GV PRO OVD is left in the eye. Special care should be taken to ensure complete removal of the Healon GV PRO OVD from the entire eye including behind the lens and the chamber angles. Continued irrigation/aspiration after displacement of the

initial bolus of viscoelastic from the eye should facilitate removal of viscoelastic if remaining in the anterior segment and behind the lens. Due to the greater viscosity of the Healon GV PRO OVD, increase in postoperative IOP may be higher than that caused by leaving the same amount of other sodium hyaluronate viscoelastic products with lower zero shear viscosity in the anterior chamber. Since rises in postoperative intraocular pressure, including cases of significant elevation and rare incidents of complications, have been reported, the following precautions are strongly recommended:

- Postoperative IOP may be elevated as a result of pre-existing glaucoma, other causes of compromised outflow, higher preoperative intraocular pressure and complications in surgical procedures (e.g. enzymatic zonulysis, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber) may also lead to increased intraocular pressure. Since the exact role of these factors is difficult to predict in any individual case, the following precautions are recommended: do not overfill the eye chambers with the Healon GV PRO OVD, completely remove the Healon GV PRO OVD by irrigation and/or aspiration at the close of the surgery and carefully monitor intraocular pressure, especially during the immediate postoperative period. Prophylactic pressure lowering treatment should always be considered, especially in cases where the Healon GV PRO OVD has to be left in the eye for clinical reasons.

- Since clinical data with use of Healon GV PRO OVD in glaucoma patients are limited, it is strongly recommended that extra care be taken to remove OVD from the eye (see below) and prophylactic pressure-lowering treatment should always be considered for these patients.
 - Start the removal directly after the IOL implantation, while the anterior chamber is still filled with the Healon GV PRO OVD and before the IOL has been centered. Go behind the IOL optic without engaging the flow of the I/A tip (port up) and then start flow. Remove the Healon GV PRO OVD from the capsular bag first and ensure that the lens has adequately centered. During removal of the Healon GV PRO OVD from the capsular bag, the continuous flow of irrigation fluid keeps the bag inflated and reduces the risk of aspirating the capsular bag. While maintaining continuous flow remove the tip from behind the optic and place it on top of the optic.
 - Continue the removal by circling the I/A tip at the iris plane, or on the optic surface, then make an additional sweep in the anterior chamber paying particular attention to the angles.

Various removal techniques may be used, depending on the surgeon's preference and/ or experience with some techniques to-remove OVDs, as advised. Descriptions for similarly difficult are provided below:

- Healon GV PRO OVD may be removed by creating maximum turbulence to make the Healon GV PRO OVD fracture into large pieces. This can be accomplished by using the described technique (described below).
 - Start by circling the hand piece in the anterior segment at iris plane.
 - Gently rest the I/A piece on the anterior surface of the optic. Press on the IOL optic on one side and rotate the I/A hand piece directing the flow into the bag. Direct the hand piece port towards the equator of the capsular bag and stay in this position for a few seconds and then repeat on the other side of the IOL optic until the Healon GV PRO OVD is completely removed. Finally, sweep the anterior chamber including the angles and repeat step 2 if necessary.

- A transient myopic shift may also occur if the Healon GV PRO OVD is not completely removed from the capsular bag behind the intraocular lens.

Using higher molecular weight, high concentration viscoelastics, such as Healon GV PRO OVD, may impede aspiration and/or potentially lead to blockage of aspiration flow through phaco tubing, especially when reuse of aspiration tubing, use of multiple vials of viscoelastics and/or lens fragments are combined.

The potential for early and short-term post-operative IOP spikes exists with cohesive OVDs, that require time and care to remove from the eye. Therefore, it is recommended that Healon GV PRO OVD be removed from the eye completely by irrigating and aspirating with sterile irrigation solution to reduce the risk of early postoperative IOP spikes.

Overfilling the anterior segment of the eye with the Healon GV PRO OVD may cause increased intraocular pressure, glaucoma, or other ocular damage.

Care should be taken to avoid trapping air bubbles behind the Healon GV PRO OVD.

The Healon GV PRO OVD is a highly purified fraction extracted from fermented bacteria which may contain minute amounts of protein

Handling instructions

Refrigerated Healon GV PRO OVD should be allowed to attain room temperature (approximately 30 minutes) prior to use.

Express a small amount of the Healon GV PRO OVD from the syringe prior to use and carefully examine it during use to avoid injecting minute rubber particles which may be released when the syringe diaphragm is punctured.

HEALON GV and HEALON are trademarks of Johnson & Johnson Surgical Vision, Inc.

Product of Sweden

Further information can be obtained from local JJV offices/sales representatives.

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