**Convenience Pack - Terminally sterilized****Bruksanvisning/DFU Healon™ PRO, CE Central**

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Ophthalmic Viscosurgical Device

### Product information

Healon™ PROviscoelastic solutionis forintraocular use. It is terminally sterilized by steam sterilization in disposable syringes delivering 0.55 or 0.85 ml respectively. The syringe is packaged together with a single use 27G cannula in a blister which is terminally sterilized by ethylene oxide.

**Composition per ml**

Sodium hyaluronate SE 10 mg
Phosphate buffer ad 1 ml

**Description**

The Healon PRO ophthalmic viscosurgical device (OVD) is a steam-sterilized, non-pyrogenic solution of a highly purified, non-inflammatory, high molecular mass fraction of sodium hyaluronate (HA) derived by microbial fermentation. The HA is provided at a concentration of 10 mg/ml (1.0%) dissolved in a physiological buffer. The Healon PRO OVD has an osmolality and a pH similar to aqueous humour. HAin theHealon PRO OVDhas an average molecular mass of 3.2 million. TheHealon PRO OVDhas a viscosity at rest about 150 Pas, 150 000 times higher than that of aqueous humour. At high shear rates, such as during injection, the viscosity of the Healon PRO OVD decreases dramatically due to shear thinning, facilitating injection through a 27G cannula. The graph represents the flow curve (shear viscosity versus shear rate).

HA is a polysaccharide made up of repeating disaccharide units linked by glycosidic bonds. HA occurs naturally on the corneal endothelium, bound to specific receptors. The natural HA layer on the corneal endothelium is reduced during irrigation, but can be restored by exogenous HA. The Healon PRO OVD has a high affinity to the receptors. The Healon PRO OVD acts as a scavenger by neutralizing the free radicals formed during cataract surgery using ultrasound.

The Healon PRO OVD is easy to inject as the viscosity decreases considerably when expelled through a cannula. TheHealon PRO OVD regards its high viscosity immediately after injection. The Healon PRO OVD creates and maintains a deep anterior chamber, allowing safe and controlled manipulation, reducing the risk of trauma to the corneal endothelium and other tissues. The Healon PRO OVD is used to atraumatically manoeuvre, control and separate tissues. The Healon PRO OVD creates and maintains a clear field of vision during surgery. The Healon PRO OVD is easy to remove from the anterior segment.

The Healon PRO OVD is eliminated in the same way as endogenous HA in the anterior segment, i.e.via the trabecular meshwork. Upon entering the general circulation HA is metabolized in the liver.

**Indication**

Eye surgery.

**Contraindications**

None known at present.

**Dosage and methods of administration**

The Healon PRO OVD is injected through a fine cannula into the eye. The amount used varies with the type of surgery, repeated injections may be required. The total amount introduced is usually less than 1 ml in anterior segment surgery.

**Undesirable effects**

A transient rise in intraocular pressure may be seen postoperatively due to a mechanical obstruction of the trabecular meshwork. The peak of this pressure rise may be expected 4-8 hours postoperatively. To minimize this risk it is recommended that the Healon PRO OVD is irrigated/aspirated out of the eye at the close of surgery. If complete removal is not achieved there is a risk of an increased intraocular pressure. Prophylactic pressure-lowering treatment should be considered in cases with glaucoma or in cases where the Healon PRO OVD for some reason is left in the eye.

**Incompatibility**

*In-vitro* studies have shown incompatibility, resulting in opalescence, between HA and solutions containing quaternary ammonium ions, e.g. detergents and benzalkonium chloride.

**Precautions**

For single use only. Do not resterilize syringe or cannula. Reuse of product or cannula may cause inflammatory reaction. Do not use if packaging has been opened or damaged. Use only the cannula provided in the carton. Do not use after date shown on product. Store cold at 2 to 8°C. Protect from light.

**Handling instructions**

Refrigerated Healon PRO OVD should be allowed to attain room temperature (approximately 30 minutes) prior to use.

Express a small amount of the Healon PRO OVD from the syringe prior to use and carefully examine it during use to avoid injecting minute rubber particles which may be released when the syringe diaphragm is punctured.

HEALON is a trademark of Johnson & Johnson Surgical Vision, Inc.

**Product of Sweden**

*Further information can be obtained from local JJV offices/sales representatives.*

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HEALON è un marchio di fabbrica di marchio di fabbrica di Johnson & Johnson Surgical Vision, Inc.