

Visual outcome, optical quality, and patients' satisfaction with a new monofocal intraocular lens, enhanced for intermediate vision: preliminary results

J Cataract Refract Surg 2020 Jan 27. Epub 2020 Jan 27. Rita Mencucci, MD, Michela Cennamo, MD, Daniel Venturi, MD, Roberto Vignapiano, MD, Eleonora Favuzza, MD

OBJECTIVES

To compare visual outcomes, contrast sensitivity, optical quality, spectacle independence and visual disturbances in patients implanted with:

- **TECNIS™** 1-piece monofocal IOL (ZCB00)
- **TECNIS Eyhance™** IOL (ICB00)

STUDY DESIGN


Prospective case series, six months follow-up. 40 patients without ocular comorbidities and corneal astigmatism <0.75 D underwent cataract surgery with bilateral implantation of the two study lenses.

- Monocular and binocular uncorrected and corrected distance (UDVA and CDVA at 4 m) visual acuities
- Uncorrected, distance-corrected and corrected intermediate (UIVA, DCIVA, CIVA at 66 cm)
- Binocular near visual acuities (UNVA, DCNVA, CNVA at 40 cm)
- Binocular defocus curve
- Photopic contrast sensitivity
- Objective scatter index (OSI), Strehl ratio, modulation transfer function (MTF) cut-off
- Halo and glare perception
- Spectacle independence


RESULTS:

TECNIS Eyhance™ IOL achieved a smoother profile along the entire curve, especially within the intermediate defocus range (up to -1.50 D, corresponding to 66 cm). The **TECNIS Eyhance™** IOL provided significantly better defocus results at -1 D and -1.5 D of defocus compared with ZCB00.

BINOCULAR DEFOCUS CURVE




SPECTACLE INDEPENDENCE FOR INTERMEDIATE DISTANCE



TECNIS Eyhance™ IOL provides higher spectacle independence for intermediate distance

CONTRAST SENSITIVITY MEASURED UNDER PHOTOPIC CONDITIONS



There is no significant difference between **TECNIS Eyhance™ IOL** and **TECNIS™ 1-piece IOL**

CONCLUSION

In patients without ocular comorbidities the **TECNIS Eyhance™ IOL (ICB00)** showed **better uncorrected intermediate visual acuity and higher spectacle independence for intermediate vision without impairment of far vision and visual quality**, compared to the **TECNIS™ IOL ZCB00**. There were **no significant differences between the two groups in terms of photopic contrast sensitivity, OSI, MTF cut-off, Strehl ratio, glare and halo perception**.

Reference:

Mencucci R, et al. Visual outcome, optical quality, and patients' satisfaction with a new monofocal intraocular lens, enhanced for intermediate vision: preliminary results. *J Caract Refract Surg* 2020 Jan 27; Epub 2020. REF2020CT4037.

©Johnson & Johnson Surgical Vision, Inc. 2021
PP2021CT6404

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.

Johnson & Johnson VISION