

# Assessment of small-bore, dual-durometer aspiration tubing in a new dual-mode Phacoemulsification System

Chetan Nirkhe M.S, Deep Mehta M.S, Jianbo Zhou PhD Johnson & Johnson Surgical Vision, Inc., Santa Ana, CA, USA

Chetan Nirkhe M.S, Deep Mehta M.S, Jianbo Zhou PhD. Johnson & Johnson Surgical Vision, Inc., Santa Ana, CA, USA


## BACKGROUND

- Bench top studies have previously shown that the aspiration tubing ID affects pump pulsation and post-occlusion surge.
- This study conducted in a lab environment used a rigid chamber with pressure sensing to replicate occlusion and measure vacuum surge levels.
- A novel aspiration tubing design was evaluated with multilayer wall thickness with variations in the material durometer to reduce compliance of the tubing while maintaining flexibility and kink resistance.

## PURPOSE

- Evaluate the performance of small-bore multilayer aspiration tubing in a new dual-pump phacoemulsification system on a benchtop setup for post-occlusion surge (POS).
- Evaluate peristaltic pump pulsation (PPP) to compare performance on different tubing designs on a benchtop setup.
- Present the venturi flow rate data using the new tubing configuration.

## POS CURVE ON PREDICATE SYSTEM TUBING DESIGNS



Reduced POS and PPP (peak to trough) with stiffer layer tubing


## POS RESULT SUMMARY-TUBING DESIGNS

Tube Type	Delta IOP -Surge (mm Hg)						Avg	% Change
	1	2	3	4	5	6		
Tube A	77.4	94.6	76.4	82.3	79.8	82.4	82.2	N/A
Tube B	71.5	64.9	70.7	71.0	71.5	71.6	70.2	14.5
Tube C	53.1	51.9	55.6	50.5	50.7	55.6	52.9	35.6

Improvements as high as 35% with multilayer stiff construction

## SMALL-BORE MULTILAYER ASP TUBING


Small diameter, two-layer PVC tubing with a soft outer thick wall and stiff inner thin wall was designed and extruded in two different configurations and used for testing on J&J dual-pump system.




## TESTING MATRIX FOR POS AND PPP

Sample Description	ID of the Asp. Tubing (mm)	Inner Layer Material	Inner Layer Thk. (mm)	Outer Layer Material	Outer Layer Thk. (mm)
Control -Single Wall Tubing (Material A)	1.21	N/A	N/A	PVC 70A (lower stiffness)	1.27
Tube Material B	1.21	PVC 90A (higher stiffness)	0.43	PVC 50A	0.84
Tube Material C	1.21	PVC 98A (highest stiffness)	0.43	PVC 50A	0.84

## POS CURVE-PREDICATE PACK& DESIGN PACK




## VENTURI MODE FLOW RATE\_20G



Reduced POS and PPP (peak to trough) on new design compared to predicate pack

## TEST SET UP DIAGRAM

Block Diagram – Fluidics Test Fixture



## CONCLUSIONS

- Using multilayer tubing design for aspiration tubing with stiff inner layer reduces POS and PPP in a benchtop study.
- With small-bore tubing, the Venturi flow rates are more controlled at higher vacuum while still achieving good balance of holdability and followability.

**Acknowledgments:** Deep Mehta, Ziad Ghannoum, Bart Johnson, Jianbo Zhou and Lelei Ji