


You have the hands of
a surgeon
and the *eyes of an artist.*

That's why your OVD is HEALON™

Trust the HEALON™ OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery¹



 **Healon™**
FAMILY of OVDs

FROM THE
MAKERS OF


TECNIS™
Family of IOLs

Healon EndoCoat™ OVD

SPEC SHEET

Ophthalmic Viscosurgical Devices (OVD)

Specifications



Healon EndoCoat™

3% Sodium hyaluronate

Healon EndoCoat™ OVD

Classification	Dispersive ¹
Composition	3% Sodium hyaluronate ¹
SKU	VT585
Volume	0.85 ml ¹
Zero Shear Viscosity	4.9 - 5.6 x 10 ⁴ cps
Average Molecular Weight	800,000 Daltons ¹
Storage Temperature	+2°C to +30°C ¹

All **HEALON™** OVDs are latex free.¹

Reference:

1. HEALON EndoCoat OVD DfU – Global multilingual – Doc. #Z311181 – Rev B – Dec 2018. REF2019CT4101.

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.