


You have the hands of  
a surgeon  
and the *eyes of an artist.*

That's why your OVD is HEALON™

Trust the HEALON™ OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery<sup>1</sup>



 **Healon™**  
FAMILY of OVDs

FROM THE  
MAKERS OF

**TECNIS™**  
Family of IOLs

# Healon GV™ PRO OVD

SPEC SHEET

Ophthalmic Viscosurgical Devices (OVD)

## Specifications



1.8% Sodium hyaluronate

### Healon GV™ PRO OVD

|                                 |                                                                                                                                                                                                                                                                                                             |
|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Classification                  | Cohesive                                                                                                                                                                                                                                                                                                    |
| Composition                     | 1.8% Sodium hyaluronate <sup>1</sup>                                                                                                                                                                                                                                                                        |
| SKUs<br>(terminally sterilized) | 10310014 – (CE Central EN, ES, PT, EL)<br>10330014 – (CE Central IT, DE, NL, FR)<br>10331014 – (CE Nordics DK, FI, SE, NO)<br>10311014 – (CE Nordics ET, LT, LV, PL)<br>10312014 – (CE East CS, HR, SR, SK)<br>10322014 – (CE East HU, RU)<br>10332014 – (CE East BG, SL, RO, TR)<br>10318014 – (CE Arabic) |
| Volume                          | 0.85 ml <sup>1</sup>                                                                                                                                                                                                                                                                                        |
| Zero Shear Viscosity            | 2.60 - 3.50 Log (Pas) <sup>2</sup>                                                                                                                                                                                                                                                                          |
| Average Molecular Weight        | 3,200,000 Daltons <sup>1</sup>                                                                                                                                                                                                                                                                              |
| Storage Temperature             | +2°C to +8°C <sup>1</sup>                                                                                                                                                                                                                                                                                   |

All **HEALON™** OVDs are latex free.<sup>1</sup>

#### References:

1. **HEALON GV PRO** OVD DfU – Doc. 52-0525-00. Rev. 2021-02-05. REF2022CT4158.
2. Evaluation of the rheological properties of **HEALON GV PRO** OVD – SEY 0389-02. REF2019CT4044.

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.