


You have the hands of
a surgeon
and the *eyes of an artist.*

That's why your OVD is HEALON™

Trust the HEALON™ OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery¹



 **Healon™**
FAMILY of OVDs

FROM THE
MAKERS OF

TECNIS™
Family of IOLs

Healon™ PRO OVD

SPEC SHEET

Ophthalmic Viscosurgical Devices (OVD)

Specifications



Healon™ PRO OVD	
Classification	Cohesive
Composition	1% Sodium hyaluronate ¹
SKUs (terminally sterilized)	10310011 (0.55 ml CE Central) 10310012 (0.85 ml CE Central) 10311012 (0.85 ml CE Nordics) 10312011 (0.55 ml CE East) 10312012 (0.85 ml CE East)
Volume	0.55 ml ¹ 0.85 ml ¹
Zero Shear Viscosity	50,000 - 562,000 mPas ²
Average Molecular Weight	3,200,000 Daltons ¹
Storage Temperature	+2°C to +8°C ¹

All **HEALON™** OVDs are latex free.¹

References:

1. **HEALON PRO** OVD DfU - CE Central - Doc. 52-0516-00 - Rev 2021-02-03. REF2022CT4156.
2. DOF2016CT0020 - Evaluation of rheological properties of **HEALON PRO** OVD and **HEALON5 PRO** OVD. 2016-06-13.

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.