


You have the hands of
a surgeon
and the *eyes of an artist.*

That's why your OVD is HEALON™

Trust the HEALON™ OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery¹



 **Healon™**
FAMILY of OVDs

FROM THE
MAKERS OF


TECNIS™
Family of IOLs

Healon5™ PRO OVD

SPEC SHEET

Ophthalmic Viscosurgical Devices (OVD)

Specifications



Healon5™ PRO

2.3% Sodium hyaluronate

Healon5™ PRO OVD

Classification	Viscoadaptive
Composition	2.3% Sodium hyaluronate ¹
SKUs (terminally sterilized)	10310015 (CE Central) 10311015 (CE Nordics) 10312015 (CE East)
Volume	0.60 ml ¹
Zero Shear Viscosity	1,300,000 - 12,600,000 mPas ²
Average Molecular Weight	3,200,000 Daltons ¹
Storage Temperature	+2°C to +8°C ¹

All **HEALON™** OVDs are latex free.¹

References:

1. **HEALON5 PRO** OVD DfU – CE Central – Doc. 52-0517-0 - Rev 2021-02-09. REF2022CT4157.
2. DOF2016CT0020 – Evaluation of rheological properties of **HEALON PRO** OVD and **HEALON5 PRO** OVD. 2016-06-13.

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.