


Sensar™ 1

1-Piece Hydrophobic Acrylic IOL

SENSAR™ 1-Piece Hydrophobic Acrylic IOL



DESCRIPTION	
OPTIC CHARACTERISTICS¹	
Powers:	+6.0 D to +30.0 D in 0.5 diopter increments
Diameter:	6.0 mm
Material:	UV-absorbing hydrophobic acrylic
Refractive Index:	1.47 at 35°C
Edge Design:	ProTEC frosted, continuous 360° posterior square edge
ULTRASOUND BIOMETRY*	
A-constant:	118.4
Theoretical AC Depth:	5.2 mm
Surgeon Factor: ²	1.45 mm
HAPTIC CHARACTERISTICS¹	
Overall Diameter:	13.0 mm
Style:	C
Material:	UV-blocking hydrophobic acrylic
Design:	Haptics offset from optic
RECOMMENDED INSERTION INSTRUMENTS	MODEL
UNFOLDER™ Platinum 1 Series Handpiece	DK7796
UNFOLDER™ Platinum 1 Series Cartridge	1MTEC30
One Series Ultra Implantation System (Syringe-style Inserter)	DK7786
One Series Ultra Cartridge	1VIPR30

*Value theoretically derived for a typical 20.0 D lens. Johnson & Johnson Vision recommends that surgeons personalize their A-constant based on their surgical techniques and equipment, experience with the lens model, and postoperative results.

To learn more and to view important safety information, please contact your Johnson & Johnson Vision Sales Representative.

References:

1. **SENSAR** 1-Piece IOL - DfU INT - Z311129, Rev. C, 09/2019. REF2021OTH4865.
2. Holladay JT. International Intraocular Lens & Implant Registry 2003. *J Cataract Refract Surg.* 2003; 29:176-197. REF2016CT0151.

©Johnson & Johnson Surgical Vision, Inc. 2021

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.

PP2021OTH6733