

TITLE: DFU, SENSAIR SOFT ACRYLIC POSTERIOR CHAMBER LENSES WITH ULTRAVIOLET ABSORBER (MULTILINGUAL, LOW COST PACKAGING)

NOTE: THE CHANGES THAT ARE MADE TO THE DFU ENGLISH PANEL CONTENT SHOULD ALSO BE APPLIED TO THE SECONDARY Z311122S LANGUAGE PANELS. DFU# Z311122P AND Z311122S NEED TO BE KEPT IN SYNC.

ISSUED

SPECIFICATIONS:

1. DIMENSIONS:

- a. Horizontal: 38-1/2" ± 1/32"
- b. Vertical: 19" ± 1/32"
- c. Flat: 19" x 38-1/2" (Tolerance: ± 1/32")
- d. Folded: 3-7/8" x 5-1/2" (Tolerance: ± 1/8")
 - i. Accordion fold 7 panels to 5-1/2" x 19"
 - ii. Right Angle fold to 5-1/2" x 11-5/8"
 - iii. Lap fold in thirds to 5-1/2" x 3-7/8"

2. STOCK: 27-30 lb. Offset or 16 lb. Bond, Smooth Opaque Finish
(Thickness range: 0.0020" – 0.0025")

3. STYLE: One sheet, printed both sides, seven panels per side, English and barcode facing out.

4. COLOR OF COPY: Black permanent pigmented ink on white paper.

5. See attached for artwork layout

6. Manufacturer specific information (i.e. Manufacture Part number) if on the label would need to be verified and confirmed by the Manufacturing site.

7. POINTS OF USE: JJSV Añasco, Puerto Rico

8. Coordinate all changes with **Z311122S**

9. DATA SHEET LAYOUT

Front	English (EN)	Graphs	Dutch (NL)	German (DE)	French (FR)	Spanish (ES)	Italian (IT)
Back	Norwegian (NO)	Danish (DA)	Finnish (FI)	Estonian (ET)	Portuguese (PT)	Swedish (SV)	Latvian (LV)

CHANGE SUMMARY

REVISION	DESCRIPTION OF CHANGE
A	Initial Release of J&J branding project.
B	Update CE Mark to MDD 93/42/EEC standards.
C	Add Ukraine conformity with notified body number. Template update to SOP501027-T020 Rev. 05

