

**TITLE: DFU, TECNIS® TORIC II 1-PIECE IOL, MODELS ZCU100-ZCU800 (INT)**

**SPECIFICATIONS:**

**ISSUED**

**1. DIMENSIONS:**

- a. Horizontal: 38-1/2" ± 1/32"
- b. Vertical: 19" ± 1/32"
- c. Flat: 19" x 38-1/2" (Tolerance: ± 1/32")
- d. Folded: 3-7/8" x 5-1/2" (Tolerance: ± 1/8")
  - i. Accordion fold 7 panels to 5-1/2" x 19"
  - ii. Right Angle fold to 5-1/2" x 11-5/8"
  - iii. Lap fold in thirds to 5-1/2" x 3-7/8"
  - iv. Two sheets folded separately.

**2. STOCK:** 27-30 lb. Offset or 16 lb. Bond, Smooth Opaque Finish.  
(Thickness range: 0.0020" – 0.0025")

**3. STYLE:** Two sheets, printed both sides, seven panels per side.

- a. Sheet 1 of 2: English faces out.
- b. Sheet 2 of 2: Hungarian faces out.

Folded sheet 2 will be placed inside of folded sheet 1 by manual assembly by the vendor.

**4. COLOR OF COPY:** Black permanent pigmented ink on white paper.

**5.** Printing to be clear and legible with no ink smears.

**6.** Insert to be clean. No visible damage, loose or attached particles (fibers) are permitted.

**7.** Vendor to package product following guidelines according to AMOS3110.

**8.** See attached for artwork layout.

**9.** Manufacturer specific information (i.e. Manufacture Part number) if on the label would need to be verified and confirmed by the Manufacturing site.

**10. POINTS OF USE:** JJSV Groningen, The Netherlands.

**11. DATA SHEET LAYOUT:**

Sheet 1: Front: English (EN), Dutch (NL), German (DE), French (FR), Spanish (ES), Italian (IT), Norwegian (NO)

Back: Danish (DA), Finnish (FI), Greek (EL), Portuguese (PT), Swedish (SV), Czech (CS), Estonian (ET)

Sheet 2: Front: Hungarian (HU), Latvian (LV), Lithuanian (LT), Polish (PL), Slovak (SK), Slovenian (SL), Serbian (SR)

Kazakh (KK)

Back: Russian (RU), Bulgarian (BG), Turkish (TR), Romanian (RO), Croatian (HR), Arabic (AR), Ukrainian (UK)

**CHANGE SUMMARY**

REVISION	DESCRIPTION OF CHANGE
XA	Initial release for Toric II lens model ZCU
A	Revision update for production release, ZCU400 was removed















