

**TITLE: DATA SHEET, FOR 1MTEC30 (MULTILINGUAL) - FOR DK7796 INSERTER BOX**

**SPECIFICATIONS:**

**ISSUED**

**1. DIMENSIONS:**

- a. Horizontal: 38.5" ± 1/32"
- b. Vertical: 24.00" ± 1/32"
- c. Flat: 38.5" x 24.00" (Tolerance: ± 1/32")
- d. Folded: 6.00" x 2.22" (Tolerance: ± 1/8")
  - i. Z fold in 12.833" x 24"
  - ii. Z fold in 6" x 12.833"
  - iii. 6" x 6.625
  - iv. 6" x 4.5"
  - v. 6" x 2.22"

**2. STOCK:** 27-30 lb. Offset or 16 lb. Bond, Smooth Opaque Finish  
(Thickness range: 0.0020" – 0.0025")

**3. STYLE:** One sheet, 2 sides, English facing out, two column text layout, double sided, 16 panels per side.

**4. COLOR OF COPY:** Black permanent pigmented ink on white paper.

**5.** Printing to be clear and legible with no ink smears.

**6.** Insert to be clean. No visible damage, loose or attached particles (fibers) are permitted.

**7.** Vendor to package product following guidelines according to AMOS #3110.

**8.** See attached for artwork layout.

**9.** Manufacturer specific information (i.e. Manufacture Part number) if on the label would need to be verified and confirmed by the Manufacturing site.

**10. POINTS OF USE:** Duckworth & Kent.

**11. DATA SHEET LAYOUT**

Side 1	English (EN)	Graphs	Dutch (NL)	German (DE)
	French (FR)	Spanish (ES)	Italian (IT)	Norwegian (NO)
	Danish (DA)	Finnish (FI)	Greek (EL)	Portuguese (PT)
	Swedish (SV)	Czech (CS)	Estonian (ET)	Hungarian (HU)

Side 2	Latvian (LV)	Graphs	Lithuanian (LT)	Polish (PL)
	Slovak (SK)	Slovene (SL)	Russian (RU)	Bulgarian (BG)
	Turkish (TR)	Romanian (RO)	Croatian (HR)	Arabic (AR)
	Serbian (SR)	Ukrainian (UK)		

**12.** Corresponding Duckworth & Kent document number is S-0328.

**CHANGE SUMMARY**

REVISION	DESCRIPTION OF CHANGE
A	Initial release for J&J branding project.
B	Add corresponding Duckworth & Kent document number. Insert barcode reference. Update lens dwell position wording and Figure 8.
C	Update CE Mark to MDD 93/42/EEC standards. Template update to SOP501027-T020 Rev. 05.
D	Update to include European Representative and to replace AMO wording







